

**HOMEFILL II
CYLINDER USAGE TIMES
(WITH PATIENT CONVENIENCE PACK)**

HomeFill II Cylinders w/PCP Conservor PULSED DOSE FLOW*					Traditional Cylinder from Delivery CONTINUOUS FLOW			
FLOW RATE	ML6	M9			M6	M9	D	E
1	9.2	14.2			2.5	3.8	5.9	10.3
2	5.0	7.6			1.3	1.9	2.9	5.1
3	3.6	5.6			0.8	1.3	1.9	3.4
4	2.9	4.5			0.6	0.9	1.5	2.6
5	2.5	3.9			0.5	0.8	1.2	2.1

* NOTE: Oxygen conservation was determined based on 20 Breaths per minute, and ideal conditions. All time estimates for HFII Patient Convenience Pack are quoted by Precision Medical based on the EasyPules O₂ Conservor.

All times quoted above are hours.