



Ameda®, Inc.

Purely Yours® Personal Double Electric Breast Pump with Tote & CustomFit™ Breast Flanges SKU # 17084



The Purely Yours Breast Pump combines **hospital-recommended technology** and **mom-friendly features**. It's perfect for mothers working full-time and for those who want a fast and easy breast-pumping experience.

Features:

Proven Airlock Protection™



- The world's only proven protective barrier.
- Helps protect breast milk and baby from bacteria, mold and viruses while pumping.¹
- Keeps tubing dry by preventing air flow between expressed milk and pump tubing while pumping.
- No need to clean the narrow tubing. Even if the container is turned, moisture and milk aerosols will not enter the tubing or pump.

CustomControl™



- Independent Speed and Suction controls make it possible for mothers to achieve a multi-phase experience.
- Each mother can customize her pump's settings to her own body's response and follow her flow to find her own best settings every time.
- Its independently adjustable controls make it possible for a mother to keep the pump set on her maximum comfortable suction level while varying the speed to more quickly trigger milk ejections

This comfortable and effective breast pump motor weighs only 1 lb. It can be conveniently powered by AC Adapter (included), 6AA batteries, or car adapter. Batteries and car adapter available separately.

Contents:

- Purely Yours electric breast pump motor unit
- Dual HygieniKit Milk Collection System includes: (2) 25.0mm breast flanges, (2) diaphragms, (4) valves, (2) tubing, (2) adapter cap, (1) white connector for single or dual pumping
- Cool 'N Carry™ Tote: (1) insulated carry bag, (6) 4 oz. bottles with lock-tight lids, (3) cooling elements
- Microfiber tote bag with adjustable shoulder strap
- CustomFit Breast Flanges, (2) 30.5mm flanges and (2) 28.5mm reducing inserts
- One-Hand Handle Assembly to convert HygieniKit to a manual pump
- AC power adapter
- Instructions, including low-literacy instructions and milk storage guidelines card
- Instructional DVD in English and Spanish

¹ Data on file

Contents on this sheet from approved Ameda marketing material 26601024-0612

Ameda, Purely Yours, HygieniKit, Proven Airlock Protection, CustomControl and associated logos are trademarks of Ameda, Inc.
© 2013 Ameda, Inc.